

Product Specifications

250V Power Supply Specifications

Input Power (Switchable)	115 VAC, 50-60 Hz 230 VAC, 50-60 Hz
Fuses	One 4A/250V, one extra fuse is provided
Output power in watts	300 watts
Output voltage range	5~250V
Output current range	10 mA~3 A
Timer	0~999.9h
Terminal pairs	4
Operating Modes	
Constant Voltage	1V step
Constant Current	1 mA step
Crossover	Auto
Display type	LCD Graphic type
Display size	53.64 x 15.64 mm (W x H)
Pause function	Yes
Safety feature	No Load Detection Load Change Detection Overload detection Ground Leak Auto Restart
Programmable	Yes Store file no.: 10 Program: up to 10 steps
Stackable	Yes
Housing material	Flame retardant ABS
Housing size	200 x 290 x 70 mm (W x D x H)
Operating temp.	0-40°C
Environmental condition	100% RH, 75 KPa-106 Kpa, Altitude not to exceed 2000 meters
Weight	2.6 kg
Certifications	CE; UL PSE VDE ROHS
Warranty	2 years