

6. Specification

Temperature Range	Ambient +5 to 100°C
Uniformity	±0.1°C
Stability	±0.2°C
Display Resolution	0.1
Heat up to 65°C -	45 minutes
Timer	Settable from 1 min. to 99 hrs. 59 min.
Orbital Radius	10mm
Linear Stroke Length	20mm
Overall dimensions (W x L x H)	33cm x 56cm x 35.6cm (13" x 22" x 14")
Weight	15.4kg (34lb)
Power Rating 115VAC - 9.65A - 60Hz.	1122w
230VAC - 4.83A - 60Hz.	1122w

7. Warranty and Service

7.1. Warranty

When used in laboratory conditions and according to these operating instructions Boekel warrants this product to be free of defective parts, material and workmanship for a period of two years from the date of shipment. The liability of Boekel Scientific for any defective equipment during the warranty period shall be limited to the repair of defective equipment or replacement thereof without charge for parts or labor.

7.2. Service

A Boekel Scientific Returned Material Authorization (RMA) number provided by Boekel Scientific is required before any Boekel products are returned for any reason. Contact Boekel Customer Service at 1-800-336-6929. A Decontamination Certificate must be completed, signed by the user, and returned to Boekel Scientific prior to receiving the RMA number. Please be sure to mark the outside of the returned goods package with this RMA number to ensure prompt handling.