
Product Description

| Models | Top Style | Top Size | Dimensions |
|----------------|------------------|-----------------|---|
| PC-200/210/220 | Ceramic | 4" x 5" | 4.4H x 5.8W x 7.3D (inches) 11.2H x 14.7W x 18.5D (cm) |
| PC-400/410/420 | Ceramic | 5" x 7" | 4.4H x 7.8W x 9.6D (inches) 11.2H x 19.8W x 24.4D (cm) |
| PC-600/610/620 | Ceramic | 10" x 10" | 4.8H x 10.5W x 14.4D (inches) 12.2H x 26.7W x 36.6D (cm) |

Operating Conditions

Corning PC-200, PC-400 and PC-600 Series Hotware is designed to be safe under the following conditions:

- Indoor use.
- Temperature of 5°C to 40°C (41°F to 104°F).
- Maximum relative humidity (non-condensing) of 80% for temperatures up to 32°C (90°F) decreasing linearity to 50% at 40°C (104°F).
- Position at least 12" from walls, 48" from ceilings and 12" between hot plates if using multiple units.
- Main supply voltage fluctuations not to exceed $\pm 10\%$ of the nominal voltage.
- Pollution Degree 2 in accordance with IEC664.

Technical Specifications - 120V

| Model | Catalog Number | Power (Volts/Hz/Watts/Amp) | Temp. Range | Stir Range | Weight |
|---------------------------|-----------------------|-----------------------------------|-----------------------|-------------------|----------------------|
| Hot Plates | | | | | |
| PC-200 | 6795-200 | 120V/60Hz/253W/2.2A | 77-1022°F 25-550°C | - | 4.0 lbs. 1.8 kg. |
| PC-400 | 6795-400 | 120V/60Hz/628W/5.3A | 77-1022°F 25-550°C | - | 6.0 lbs. 2.7 kg. |
| PC-600 | 6795-600 | 120V/60Hz/1043W/8.7A | 77-1022°F 25-550°C | - | 10.0 lbs. 4.5 kg. |
| Stirrers | | | | | |
| PC-210 | 6795-210 | 120V/60Hz/33W/0.3A | - | 60-1100 | 4.5 lbs. 2.0 kg. |
| PC-410 | 6795-410 | 120V/60Hz/73W/0.7A | - | 60-1100 | 6.5 lbs. 2.9 kg. |
| PC-610 | 6795-610 | 120V/60Hz/73W/0.7A | - | 60-1100 | 10.0 lbs. 4.5 kg. |
| Stirrer/Hot Plates | | | | | |
| PC-220 | 6795-220 | 120V/60Hz/283W/2.4A | 77-1022°F 25-550°C | 60-1100 | 5.0 lbs. 2.3 kg. |
| PC-420 | 6795-420 | 120V/60Hz/698W/5.9A | 77-1022°F 25-550°C | 60-1100 | 7.0 lbs. 3.2 kg. |
| PC-620 | 6795-620 | 120V/60Hz/1113W/9.3A | 77-1022°F 25-550°C | 60-1100 | 11.5 lbs. 5.2 kg. |