

## Section 2 Specifications

Operating Range .....30° to 75°C

**Note** Operating temperature must be at least 5°C above ambient. ▲

Average Uniformity @ 37°C\* .....±0.3°C

Resolution .....0.1°C

Control Sensitivity .....±0.1°C

Recovery Time @ 37°C (door opened 30 seconds)

Small incubator .....2.0 minutes

Medium incubator .....2.0 minutes

Large incubator .....3.0 minutes

*\*When tested per ASTM E-1292-94*

### Operating Environment

Temperature Range .....10° to 35°C

Humidity Range (non-condensing) . . 0% to 90% RH

### Electrical Requirements (includes 5 amp outlet)

#### Small incubators

11-690-625F .....120V, 60 Hz, 8.3A, 996W

11-690-626F ... .240V, 50/60 Hz, 6.6A, 1560W

#### Medium incubators

11-690-637F .....120V, 60 Hz, 8.3A, 996W

11-690-638F ... .240V, 50/60 Hz, 6.16A, 1560W

#### Large incubators

11-690-650F .....120V, 60 Hz, 8.9A, 1068W

11-690-651F ... .240V, 50/60 Hz, 7.0A, 1680W

**Chamber Volumes**

Small incubators . . . . .	2.2 cu ft
Medium incubators . . . . .	3.4 cu ft
Large incubators . . . . .	4.5 cu ft

**Chamber Dimensions (W x D x H)**

Small incubators . . . . .	18 x 18 x 12.5 in
Medium incubators . . . . .	18 x 18 x 19 in
Large incubators . . . . .	18 x 18 x 25.5 in

**Environmental Conditions**

Indoor use only

Operating: 17°C to 27°C; 20% to 80% RH, non-condensing.

Installation Category II (overvoltage) in accordance with IEC 664.

Pollution Degree 2 in accordance with IEC 664.

Altitude Limit: 2,000 meters.

Storage: -25°C to 65°C, 10% to 85% relative humidity (RH)