

# Appendix B

## Product Specifications

<b>iBind™ Flex Western Device specifications</b>	Dimensions:	30.0 cm (l) × 25.2 cm (w) × 8.0 cm (h)
	Material:	Aluminum, plastic (PC/ABS), steel, silicone, neodymium (magnets)
	Operating Temperature:	18°C to 30°C
	Temperature Limit:	30°C

The iBind™ Flex Western Device is impervious to alcohol, but not compatible with chlorinated hydrocarbons (e.g., chloroform), aromatic hydrocarbons (e.g., toluene, benzene) or acetone.

<b>iBind™ Flex Card specifications</b>	Dimensions:	17.8 cm (l) × 17.8 cm (w) × 0.8 cm (stack height)
	Material:	Glass fiber
	Operating Temperature:	18°C to 30°C
	Temperature Limit:	30°C