



Overview

Product Shelves accommodate bulk samples, microtiter plates, vials or serum bottles. They are designed to fit inside Drying Chambers (sold separately).

Three removable, 10.0" diameter aluminum shelves with 1/4" rim. Shelves slide into stainless steel support. Includes three rubber feet. Removable shelf design allows for flexible shelf height spacing. Shelves can be placed in freezer to pre-freeze sample for quick loading.

Catalog Number: 7317100

Specifications

- **Weight:** 4.5 lbs
- **Weight metric:** 2.5 kg
- **Dimensions:** 10.9" w x 10.4" d x 9.8" h
- **Dimensions metric:** 27.7 x 26.4 x 24.9 cm
- **Applicable Products:** 16-Port Drying Chamber, Large Clear Chamber, Large Clear Chamber with Valves, Short Clear Chamber, Short Clear Chamber with Valves, Tall Clear Chamber, Tall Clear Chamber with Valves

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

