

Acid Resistant High Capacity Micro Rotor



View online: <https://www.labconco.com/product/high-capacity-micro-rotor-ptfe-coated/5787>



Overview

Catalog Number: 7553501

Specifications

- **Weight:** 2.5 lbs
- **Weight metric:** 1.1 kg
- **Dimensions:** 12"x12"x12"
- **Applicable Products:** CentriVap Concentrators
- **Chemical Resistance:** Acid-Resistant
- **Material:** PTFE-coated aluminum

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

