

F14-6x250LE Fixed Angle Carbon Fiber Rotor

Rotor Specifications

Capacity (ml)	6 x 250
Tube Dimensions (mm)	62 x 125
Angle °	23
Maximum Speed (rpm)	11,000*
K Factor	2578
Net Weight (kg)	6,6

* Max Speed may vary depending on centrifuge unit.

RCF/Radius	RCF (x g)	Radius (cm)
Maximum	17.180	12,7
Minimum	5.001	3,7



F14-6x250LE

Cat. No.	Description
75003662	F14-6x250LE Carbon Fiber Rotor

Accessories

Cat. No.	Description
099-062153	Lid Assembly
021-062153	O-ring Refresher Kit

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims ØxL (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:		
							Cat. No.	Qty /Pkg	Description	Cat. No.	Qty /Pkg	No. of Places	Cat. No.	Qty /Pkg	No. of Places
3141-0250	250	250	PP Oak Ridge Bottle	4	11.000	61x134	Included	4	PP Sealing	-	-	-	-	-	-
3140-0250	250	250	PC OakRidge Bottle	4	11.000	61x134	Included	4	PP Sealing	-	-	-	-	-	-
334959	50	50	Nunc Conical	25	-	30x121	Included	25	PP Sealing	03072	1	1	-	-	-
3139-0050	50	50	PP Oak Ridge Tube	10	11.000	29x114	Included	10	PP Sealing	010-0138	2	1	-	-	-
3138-0050	50	50	PC Oak Ridge Tube	10	11.000	29x114	Included	10	PP Sealing	010-0138	2	1	-	-	-
3114-0050	50	50	Teflon Oak Ridge Tube	2	11.000	29x114	Included	2	PP Sealing	010-0138	2	1	-	-	-
-	30	-	Glass Tube	-	-	24x105	-	-	-	00449	1	3	-	-	-
366036	15	15	Nunc Conical	50	-	17x120	Included	50	PP Sealing	00456	1	5	-	-	-
-	10	-	Blood Collection	-	-	16X100	-	-	-	00389	1	7	-	-	-
3139-0010	10	8	PP Oak Ridge Tube	10	11.000	16x82	Included	10	PP Sealing	00389	1	7	-	-	-
3138-0010	10	8	PC Oak Ridge Tube	10	11.000	16x82	Included	10	PP Sealing	00389	1	7	-	-	-
-	7	-	Blood Collection	-	-	13x100	-	-	-	00369 + 00456	1	5	-	-	-
-	5	-	Blood Collection	-	-	13x75	-	-	-	00388	1	12	-	-	-
-	5	-	RIA Tubes	-	-	13x75	-	-	-	00388	1	12	-	-	-