

ELECTRONIC BENCHTOP AUTOCLAVES

Designed to autoclave wrapped, unwrapped, glass and tools. Not suitable for liquids (see liquid autoclaves on pages 8 - 12).

- Sterilization time up to 99 minutes
- Temperature selection from 105° to 137°C (221° to 279°F)
- Includes trays
- EAP Models allow closed door drying
- In-line particle filter to protect valves
- Electronic pressure and temperature sensors
- Low Water Alarm - sounds when reservoir needs to be refilled.
- Standby Heating - keeps chamber warm for quicker cycles
- Built In Printer - Provides a printed summary at the end of each cycle, documenting temperature, elapsed sterilization and drying time, selected cycle number, real time and date, as well as any errors or interruptions that may occur during the cycle.



	2340EP	2540EAP
Chamber Volume (L)	19	23
Max Available Height (in)	6.75"	8.3"
Chamber Dimension (ØxD)	9" x 18"	10" x 18"
External Dimension (WxHxD)	20" x 15" x 20"	20" x 15" x 20"
Power	120V	120V
Utilities	No Plumbing Required	
P/N	023210215	023210380

	3545EP	3850EP	3870EAP
Chamber Volume (L)	34.4	65	85
Max Available Height (in)	10.1"	12.1"	12.1"
Chamber Dimension (ØxD)	12.2" x 18.5"	15" x 20"	15" x 27"
External Dimension (WxHxD)	23.5" x 18" x 21.5"	26" x 21" x 27"	26" x 21" x 34"
Power		230V, 1 phase (60 Hz) NEMA 6-15	
Utilities	No Plumbing Required		
P/N	023210690	023210495	023210673
Autoclave Stand (PN 023212609)	NA	Free with Autoclave	Free with Autoclave
Table with SS top (30" x 36") & Castors	023212601	023212601	023212601

2 YEAR WARRANTY