

PARTS LIST 1326 & 1327F

Description	115V	220V
Main Temperature Control	1750863	1750863
Over Temperature Control	1750615	1750648
Cord Set	1800516	1800537
Door Handle	3800610	3800610
Heating Element 1326	9570774	9570803
Heating Element 1327F	9570778	N/A
On/Off Switch	7850570	7850570
Knob Main Temperature	4450528	4450528
Knob Over Temperature	4450506	4450506
Pilot Light Green	4650554	4650554
Pilot Light Red	4650553	4650553
Shelf Clips	1250511	1250511
Shelf	5130714	5130714
Thermometer Clip	5080865	5080865
Thermometer	8200509	8200509
Fan Motor 1325F	4880527	4880528
Fan Blade 3" 1325F	2600545	2600545
Fan Blade 4.75" 1325F	2600502	2600502
Filter EMI	2800502	2800502
Fuse	3300516	3300515
Inlet with Fuse Drawer	4200505	4200505
Fuse Holder	N/A	3300501

UNIT SPECIFICATIONS

Weight	Shipping	Net
1326	153 lbs.	113 lbs.
1327F	153 lbs.	113 lbs.

Dimensions	Exterior WxDxH	Interior WxDxH
1326	24 X 24 X 37"	20 X 19 X 28"
1327F	24 X 24.5 X 37"	20 X 17.5 X 28"

Capacity	Cubic Feet
1326	6.1
1327F	5.7

Temperature	Range	Uniformity
1326	Amb+5° to 200°C	<u>±</u> 4°C
1327F	Amb+5° to 200°C	<u>±</u> 4°C