

SPECIFICATIONS

MODEL 1000

Overall dimensions LxWxH		11 x 17 x 4" (27.9 x 43.2 x 10.2cm)
Tray dimensions LxW		8.75 x 11.75" (22.2 x 29.8cm)
Electrical	120 volts:	5 amps, 25 watts
	230 volts:	2.5 amps, 25 watts
Fuses		5mm x 20mm, 5 amp quick acting
Speed Range		40 to 300 rpm
Speed Accuracy	Advanced:	
Above 100rpm		±2% of set speed
Below 100rpm		±2rpm
Orbit		0.6" (15 mm)
Capacity		8 lbs (3.6kg)
Controls	Standard:	See page 6
	Advanced:	See page 8
Ship Weight		25 lbs (11.3kg)

MODEL 3500

Overall dimensions LxWxH		16.25 x 14 x 5.75" (41.3 x 35.6 x 14.6cm)
Tray dimensions LxW		13 x 11" (33 x 28cm)
Electrical	120 volts:	5 amps, 75 watts
	230 volts:	2.5 amps, 75 watts
Fuses		5mm x 20mm, 5 amp quick acting, 250V
Speed Range	Advanced:	15 to 500 rpm
	Standard:	25 to 500 rpm
Speed Accuracy	Advanced:	
Above 100rpm		±1% of set speed
Below 100rpm		±1rpm
Orbit		0.75" (19mm)
Capacity		~35 lbs (16kg) @ 75 rpm ~5 lbs (2.3kg) @ 500 rpm
Timer	Advanced:	1 second to 160 hours
	Standard:	1 to 120 minutes
Controls	Advanced:	See page 8
	Standard:	See page 6
Ship Weight		49 lbs (22.2 kg)

MODEL 3750

Overall dimensions LxWxH		16.25 x 14 x 5.75" (41.3 x 35.6 x 14.6cm)
Tray dimensions LxW		13 x 11" (33 x 28cm)
Electrical	120 volts:	5 amps, 40 watts
	230 volts:	2.5 amps, 40 watts
Fuses		5mm x 20mm, 5 amp quick acting, 250V
Speed Range		20 to 300 rpm
Speed Accuracy		
Above 100rpm		±2% of set speed
Below 100rpm		±2rpm
Orbit		0.75" (19mm) Total Travel
Capacity		~15 lbs @ 300 rpm
Timer		1 second to 160 hours
Controls		See page 8
Ship Weight		49 lbs (22.2 kg)

MODEL 5000

Overall dimensions LxWxH		23.6 x 26.7 x 6.0" (59.9 x 67.8 x 15.2cm)
Tray dimensions LxW		18 x 24" (45.7 x 61cm)
Electrical	120 volts:	5 amps, 75 watts
	230 volts:	2.5 amps, 75 watts
Fuses		5mm x 20mm, 5 amp quick acting, 250V
Speed Range	Advanced:	20 to 500 rpm
	Standard:	25 to 500 rpm
Speed Accuracy	Advanced:	
Above 100rpm		±1% of set speed
Below 100rpm		±1rpm
Orbit		1" (25.4 mm)
Capacity		~50 lbs @ 125 rpm ~10 lbs @ 500 rpm
Timer	Advanced:	1 second to 160 hours
	Standard:	1 to 120 minutes
Controls	Advanced:	See page 8
	Standard:	See page 6
Ship Weight		109 lbs (49.4 kg)