

UNIT SPECIFICATIONS

Weight	Shipping	Net
5015	195 lbs.	134 lbs.
5025	400 lbs.	268 lbs.

Dimensions	Exterior WxDxH (in.)	Interior WxDxH (in.)
5015	25.75 x 28.5 x 37	18 x 20 x 25.25
5025	25.75 x 28.5 x 73.5	*18 x 20 x 25.25

Capacity	Cubic Feet
5015	6
5025	*6

Temperature	Range	Control	Uniformity	CO2 Range
5015	5° above amb. to 70°C	±0.1°C	±0.25° @ 37°C	0-20%
5025	5° above amb. to 70°C	±0.1°C	±0.25° @ 37°C	0-20%

* each chamber