

SPECIFICATIONS - ANALOG MULTI-TUBE VORTEXER

Overall dimensions (L x W x H):	9.5 x 15.1 x 16" (24.1 x 38.4 x 40.6cm)
Electrical (50/60 Hz):	100 - 120V - 0.8 amps, 100 watts 220 - 240V - 0.4 amps, 100 watts
Fuses:	120V - 5mm x 20mm, 5 amp, 250V 230V - 5mm x 20mm, 5 amp, 250V
Fuses:	quick acting
Speed range:	1200 to 2400rpm
Accuracy:	+/- 25rpm
Orbit:	3.6mm
Weight capacity:	10lbs (4.5kg)
Timer:	mechanical timer runs 0 to 60 seconds
Controls:	rocker switch: run, standby, timer, speed knob, variable 1 to 10 dial markings, time knob, variable 1 to 10 dial markings
Ship weight:	45lbs (20.4kg)

