

Product Name:	HCUCFS-0420
Description:	<p>VWR Catalog No. 10819-894 VWR Undercounter Freestanding Freezer 4 Cu. Ft. Capacity One solid swing door, right hinged Mechanical thermostat Adjustable operating temperature range: -15°C – -25°C 2 fixed shelves Keyed door lock Magnetic door gasket for positive seal Recessed door handle Leveling legs Probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 20 1/4"W x 21 5/8"D x 31 1/2"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, .90 Amps, 1/7 HP UL/C-UL listed Shipping weight: 94 lbs.</p>