

The Waters® ACQUITY UPLC™ FlexCart wheeled cart is a complete system platform that facilitates convenient installation of the ACQUITY UPLC and makes it easier to position the system close to the inlet of the mass spectrometer. The FlexCart provides electrical outlets, supports computer monitor and keypad, and manages fluid waste container.

The ACQUITY UPLC FlexCart has maximum accessibility from the sides for workstation and waste-can interface. AC power is available at the rear of the cart for distribution to the requisite ACQUITY UPLC modules and peripherals.

The FlexCart's mobile platform allows the ACQUITY UPLC system to be moved among several workplaces. Its adjustable tabletop height allows optimum positioning of the column outlet of the system with a mass spectrometer.



Specifications

Dimensions deep	30 inches (76.1 cm) wide x 33 inches (83.8 cm)
Minimum Table Height	29.1 inches (73.9 cm)
Maximum Table Height	41.5 inches (105.4 cm)
Height Adjustment	12.3 inches (31.2 cm)
Cart Weight	170 lbs (77.3 kg)
Supportable Weight	400 lbs (181.8 kg)

Order Information

ACQUITY UPLC FlexCart

Part Number

205015015

Waters

For Complete  Confidence